

## General Surgery

### Taper Point Needles

Ethicon	USS	D&G	Length
<b>1/2 CIRCLE</b>			
FN-2, CT-3	GS-23	T-3	22 mm
FN, CT-2	GS-22	T-19	27 mm
CT-1	GS-21	T-12	37 mm
CT	GS-24	T-25	40 mm
CTX	GS-25	T-20	48 mm
TP-1, XLH <sup>†</sup>	GS-26	T-60	65 mm
XLH <sup>†</sup>	GS-27	T-29	76 mm
<b>5/8 CIRCLE</b>			
UR-6	GU-46	TT-2, TT-20	27 mm
UR-5	GU-45	TT-3	37 mm
UR-4	GU-44	TT-4	40 mm
<b>1/2 CIRCLE</b>			
MO-7	HGS-23	T-3	22 mm
MO-6	HGS-22	T-40	27 mm
MO-5	HGS-20	T-41	34 mm
MO-4	HGS-21	T-56	37 mm
MO-2	HGS-24	T-25 <sup>†</sup>	40 mm

### Blunt Point Needles

Ethicon	USS	D&G	Length
<b>1/2 CIRCLE</b>			
CTB-1	BGS-21	MT-12	37 mm
-	BGS-24	MT-25	40 mm
CTX-8	BGS-25	MT-20	48 mm
BP-1	BGS-29	MT-60 <sup>†</sup>	63 mm
-	BGS-27	BT-60	64 mm
BP	BGS-28	-	85 mm
<b>3/8 CIRCLE</b>			
-	-	NE-9	64 mm

### Straight Taper Point Needles

Ethicon	USS	D&G	Length
-	-	TS-4	44 mm
ST	ST-1	-	51 mm
ST-1	ST	TS-9	60 mm

## Cuticular

### Reverse Cutting Needles

Ethicon	USS	D&G	Length
<b>3/8 CIRCLE</b>			
C-2	C-1	CE-2 <sup>†</sup>	11 mm
FS-3	C-12	CE-3	16 mm
FS-2	C-13	CE-4	19 mm
FS-1	C-14	CE-6	24 mm
FS	C-15	-	26 mm
FSL	C-16	CE-10	30 mm
FSLX	C-17	CE-12	39 mm
LR	GS-18	CE-24	77 mm
<b>1/2 CIRCLE</b>			
M-1	C-21	-	13 mm
-	C-22	C-4 <sup>†</sup>	20 mm
X-1	C-23	C-30 <sup>†</sup>	24 mm
OS-4	C-25	-	30 mm

### Straight Cutting Needles

Ethicon	USS	D&G	Length
-	SC-11	-	11 mm
-	SC-4	-	24 mm
-	SC	CS-3	44 mm
TS	SC-1	-	51 mm
KS	SC-2	CS-1	60 mm
-	SC-10	CS-2	81 mm
KS-3, SKS-3	SC-6	CS-6	88 mm
-	-	CS-10	250 mm

Not Actual Size

## Microsurgery

### Taper Point Needles

Ethicon	USS	D&G	Length
<b>3/8 CIRCLE</b>			
BV-50-3/2 <sup>†</sup>	MV-50-3	TE-50	4 mm
BV-75-3 <sup>†</sup>	MV-70-3	TE-70	4 mm
BV-100-3	MV-100-3	TE-143	4 mm
BV-130-3	MV-135-3	TE-143	4 mm
BV-75-4	MV-70-4	TE-100	5 mm
BV-100-4	MV-100-4	TE-100	5 mm
BV-130-4/5 <sup>†</sup>	MV-135-4	CV-345, TE-100	5 mm
BV-130-5	MV-135-5	CV-345, TE-145	6 mm
BV-175-6	MV-175-8	CV-351	8 mm
BV-175-6	MVF-175-8	CV-351	8 mm
BV-175-8	MV-175-9	TE-11, TE-10	10 mm

## Cardiovascular

### Taper Point Needles

Cardiovascular, Gastrointestinal, Pediatric

Ethicon	USS	D&G	Length
<b>3/8 CIRCLE</b>			
BV-1	CV-1, CVF-1	TE-10 <sup>†</sup> , TE-11 <sup>†</sup>	9 mm
BV	CV, CVF	TE-9	11 mm
C-1	CV-11, CVF-11	TE-1 <sup>†</sup>	13 mm
BB	CV-15	TE-7 <sup>†</sup>	17 mm
BB-1	CV-13	-	22 mm
TE	CV-17	TE-3 <sup>†</sup>	32 mm
-	CV-19	-	50 mm
<b>1/2 CIRCLE</b>			
RB-3 <sup>†</sup>	CV-20	T-37, CV-337	10 mm
TF	CVF-21	T-30 <sup>†</sup>	12 mm
RB-2	CV-22, CVF-22	T-30	13 mm
RB-1	CV-23, CVF-23	T-31, CV-331 <sup>†</sup>	17 mm
SH-2	CV-24	T-16	20 mm
SH-1	CV-25	-	22 mm
SH	V-20, VF-20	T-5, CV-305 <sup>†</sup>	26 mm
-	V-30	T-15	30 mm
MH	V-26	T-10	37 mm

### Penetrating Point (DT Series)

### Tapercutting Needles (KV Series)

Vascular & Pediatric

Ethicon	USS	D&G	Length
<b>3/8 CIRCLE</b>			
CC	KV-1, KV-1	DTE-10	9 mm
CC-1	KV-11, KV-11	DTE-1, DTE-2 <sup>†</sup>	13 mm
V-4	KV-15	DTE-7 <sup>†</sup>	17 mm
<b>1/2 CIRCLE</b>			
V-5	KV-5,	DT-31	17 mm
V-7	KV-7	DT-5	26 mm
<b>General Closure-Sternum</b>			
<b>1/2 CIRCLE</b>			
V-34	KV-34	DT-12	37 mm
V-37	KV-37	DT-25	40 mm
V-40	KV-40	DT-20	48 mm
V-56 <sup>†</sup>	KV-56	-	60 mm

### Tapercutting Needles

Ethicon	USS	D&G	Length
<b>3/8 CIRCLE</b>			
-	MVK-70-3	CTE-103	4 mm
V-100-3, V-130-3 <sup>†</sup>	MVK-100-4	CTE-100	5 mm

## Ophthalmic

### Premium Point Spatula Needles

Ethicon	USS	D&G	Length
<b>1/4 CIRCLE</b>			
CS/TG-100-8	SE-100-8	LO-90-8	6 mm
<b>3/8 CIRCLE</b>			
CS/TG-140-6	SE-140-6	LE-140-6 (LE-1)	6 mm
CS/TG-140-8	SE-140-8	LE-140-8 (LE-2)	6 mm
<b>1/2 CIRCLE</b>			
CS/TG-160-4	SE-160-4	LA-160-4 (LA-100)	6 mm
CS/TG-160-6	SE-160-6	LA-160-6 (LA-1)	6 mm
TGW-160-6	-	-	6 mm
CS/TG-160-8	SE-160-8	LA-160-8 (LA-2)	6 mm
CS/TG-175-8	SE-175-8	L-175-8 (L-2)	7 mm
<b>COMPOUND CURVE</b>			
CS-C-6 <sup>†</sup>	SE-CC-6	LB-90/50-6	5 mm

### Spatula Needles

Ethicon	USS	D&G	Length
<b>1/4 CIRCLE</b>			
S-29	SS-29	DO-3	7 mm
S-24	SS-24	DO-5, DO-1	8 mm
S-14	SS-14	DO-2	10 mm
<b>1/2 CIRCLE</b>			
S-2	SS-2	D-1	8 mm

### Premium Point Reverse Cutting Needles

Ethicon	USS	D&G	Length
<b>3/8 CIRCLE</b>			
G-6	HE-6	CE-20	7 mm
G-1	HE-1	CE-21	11 mm
G-3	HE-3	CE-23	13 mm
<b>1/2 CIRCLE</b>			
G-7	HE-7	C-21	7 mm
G-2	HE-2	C-22	9 mm

## Endoscopic

### Taper Point Needles

Ethicon	USS	D&G	Length
-	ESK	TH-2	23 mm
ST-4	EST	TS-20	19 mm

(Cardiovascular Needles continued)

### Conventional Cutting Needles

Ethicon	USS	D&G	Length
<b>1/2 CIRCLE</b>			
-	-	SC-4	21 mm
CCS	SCC	CV-320	48 mm
CCS-1	SCC-1	-	60 mm
<b>3/8 CIRCLE</b>			
-	-	GCC-90 COLTS	90 mm

## Plastic/Cosmetic

### Premium Reverse Cutting Needles

Ethicon	USS	D&G	Length
<b>3/8 CIRCLE</b>			
P-6	P-16	PRE-20	7 mm
P-1	P-10	PRE-1, SBE-1	11 mm
P-3	P-13	PRE-2, SBE-2, PRE-302	13 mm
PS-3	P-11	PRE-3, SBE-3	16 mm
PS-2	P-12	PRE-4, SBE-4, PRE-304	19 mm
PS-1	P-14	PRE-6, SBE-6	24 mm
<b>1/2 CIRCLE</b>			
P-2	P-21	PR-1, C-1	9 mm
PS-5	P-22	PR-2	13 mm
PS-4	P-24	PR-4 <sup>†</sup>	16 mm

### Premium Conventional Cutting Needles

Ethicon	USS	D&G	Length
<b>3/8 CIRCLE</b>			
PC-1 <sup>†</sup>	PC-13	PE-1 <sup>†</sup>	11 mm
PC-3	PC-11	PE-3	16 mm
PC-5	PC-12	PE-5	19 mm

## Orthopedic

### Reverse Cutting Needles

General Closure

Ethicon	USS	D&G	Length
<b>1/2 CIRCLE</b>			
CP-2	GS-10	C-6	26 mm
CP-1	GS-11	C-9	37 mm
CP	GS-12	C-23	40 mm
CPX	GS-13	C-12	48 mm
<b>Orthopedic, OB/GYN</b>			
<b>1/2 CIRCLE</b>			
OS-4	HOS-10	C-40	26 mm
OS-6	HOS-11	C-41	37 mm
OS-8	HOS-12	C-42	40 mm

(Endoscopic Taper Point Needles continued)

Ethicon	USS	D&G	Length
<b>1/2 CIRCLE</b>			
SH-2	EV-23	T-31	17 mm
SH-1	EV-25	T-16	22 mm
SH	EV-20	T-5, CV-305	26 mm
CT-2	EGS-22	T-19	27 mm
<b>1/2 CIRCLE</b>			
ST-1	EKS	-	24 mm

### Straight Tapercutting Needles

USS D&G Product	Material	Ethicon
POLYSORB*	Coated, Braided Synthetic Absorbable	Vicryl (Coated)
DEXON* II DEXON* S	Coated and Uncoated Braided, Synthetic Absorbable	Vicryl (Coated)
BIOSYN*	Monofilament, Synthetic Absorbable	Microneedle
MAXON*	Monofilament, Synthetic Absorbable	PDS II
MILD & CHROMIC GUT	Natural Absorbable Gut	Surgical Gut
PLAIN GUT	Natural Absorbable Gut	Surgical Gut
SURGIPO* SURGILENE*	Nonabsorbable, Monofilament Polypropylene	Prolene
NOVAFIL*	Nonabsorbable, Monofilament Polybutester	Prolene
VASCUFIL*	Nonabsorbable Monofilament Polybutester with Polytribolate Absorbable Coating	Pronova
SURGIDAC* TI-CRON*	Nonabsorbable Coated Braided Polyester	Mersilene (Uncoated) Ethibond (Coated)
MONOSOF* DERMALON*	Nonabsorbable Monofilament Nylon	Ethilon
BRALON* SURGILON*	Nonabsorbable Coated Braided Nylon	Nylon
SOFSILK* SILK	Nonabsorbable Coated Braided Silk	"Permahand" Silk (Coated)
STEEL	Nonabsorbable 316L Stainless Steel	Surgical Stainless Steel
FLEXON*	Temporary Cardiac Pacing Wire	Surgical Stainless Steel

tyco  
Healthcare

150 Glover Avenue, Norwalk, CT 06856 USA  
1 800 722 USSC · FAX: 1 800 544 USSC  
©2001 United States Surgical, a division of Tyco Healthcare Group LP  
All Rights Reserved. Printed in USA, PP 7.5M 2/01 S50000011  
<sup>†</sup>Suggested Replacement

United States  
Surgical